


# MW-12A

## Microplate Washer

### Technical Specifications

Wash programs	up to 100
Wash cycles	0 to 9
Dispensing volume	50 up to 3000 µL
Bottom collision protection	automatic
Residual volume per well	< 2 µL
Soaking time	0 to 300 seconds
Shaking time	0 to 300 seconds
Multiple options	1 x 8, 1 x 12
Standard bottle set	2 x 2L (detergent), 1 x 2L (Waste) 1 x 2L (distilled water)
Power supply	AC 100-240V, 50/60Hz
Temperature	5°C~35°C
Humidity	15~85%
Dimension	400mm (L) x 330mm (W) x 145mm (H)
Weight	7.2 Kg



# MW-12A

## Microplate Washer

- Compatible with flat, U- or V-bottom plates
- 2 wash modes for 8- and 12-well strips
- Automatic bottom detection and collision protection
- Variable Shaking and soaking time
- Up to 100 user-defined protocols
- Sweep function with aspiration
- On-board automatic self-diagnostic program
- Automatic liquid monitoring system



Mindray Building, Keji 12th Road South,  
High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R. China  
Tel: +86 755 8188 8998 Fax: +86 755 26582680  
E-mail: intl-market@mindray.com www.mindray.com

Mindray is listed on the NYSE under the symbol "MR"

mindray is a trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd.  
© 2013 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved. Specifications subject to changes without prior notice.  
P/N: ENG-MW12A-210285x2-20130303



**mindray**

**mindray**  
healthcare within reach